



Spectra Iris

Your travel-sized BIO

Spectra Iris

From making rounds at your hospital to serving remote communities, Spectra Iris was designed to go wherever your work takes you. Lightweight and portable, Spectra Iris offers reputable Keeler Optics and brilliant LED illumination for a crystal-clear indirect exam.

Key features

- **Keeler optics** - we incorporate high-quality optics in all binocular indirects, and Spectra Iris is no exception. Precision-crafted lenses and mirrors provide an unhindered stereoscopic fundus view.
- **Spectacle-mounted** - Spectra Iris attaches to a stylish sports frame, allowing you to take it on and off with ease. And best of all, no more helmet hair!
- **Lightweight & portable:** by condensing the optical system and controls, Spectra Iris offers unmatched convenience. Slip it into your pocket or wear it around your neck wherever you go.



New carrying case

Sling it over your shoulder or carry it around - Spectra Iris now includes a custom-designed, portable carrying case. A perfect accessory for protecting your indirect and condensing lens on-the-go.*

**Condensing lens is an optional accessory.*

Part numbers

Spectra Iris kit

1205-P-2010 Spectra Iris kit; small frame

1205-P-2011 Spectra Iris kit; medium frame

1205-P-2012 Spectra Iris kit; large frame

Optional accessories

3412-P-5354 Spectra Iris carrying case

1205-P-7000 Spectra Iris teaching mirror

1919-P-5282 Spectra Iris battery pack

1941-P-5414 Spectra Iris lithium charger

2105-K-1159 Black Volk 20D lens

What's in the box?

Spectra Iris kit

- Spectra Iris indirect ophthalmoscope Keeler
- black sport frame (S/M/L depending on part number)
- Spectra Iris carrying case
- Spectra Iris lithium charger
- Transformer with plug adapters prescription
- lens holder
- Lens cloth
- Lanyard
- Screwdriver
- Spectra Iris instructions for use

Specifications

Spectra Iris

PD range	48-76 mm
Filters	<input type="radio"/> Neutral <input checked="" type="radio"/> Cobalt blue <input type="radio"/> Red-free
Spot size	20 mm to 60 mm continuous
Dimensions	7.5 x 0.5 x 12 cm (H x W x D)
Weight	136g (with sports frame)
Illumination	LED
Colour Temperature	4000K

Battery

Type	Lithium-ion
Duration	4 hours continuous at maximum brightness
Weight	170g

Carrying case

Dimensions	31 x 22 x 12 cm (H x W x D)
-------------------	-----------------------------



More Keeler instruments

With over 100 years of industry experience, Keeler has become a household name among optometrists and ophthalmologists worldwide.

View our full range at

www.keeler.co.uk | www.keelerusa.com

Got a question? - call us

UK: +44 (0)1753 857177 India: +91 99 3031 1090

USA: +1 800 523 5620 Brazil: +55 11 4302 6053

China: +86-18512119109



A Halma company

Manufacturing address



Keeler Ltd. Clewer Hill Road Berkshire SL4 4AA UK



053122 REV D