# PSL classic

Excellence in optics





New from Keeler. Please see reverse for details.





## Keeler – British manufacturing experience and quality delivering great value

#### Precision machined aluminium chassis

Keeler designed and manufactured with the strength to withstand the daily traumas of the demanding ophthalmic environment. Lightweight, durable and robust.

#### Advanced optics, x10 & x16 magnifications

Top quality multi layered glass is crafted into our Keeler multi-element and multi-layer coated optical system to give maximum performance. Clarity is guaranteed.

#### Controllable illumination

Proven halogen illumination with illumination levels controlled by a rheostat, giving you total control of your light levels. From zero to maximum, you choose.

#### Big Slit Lamp features, Portable usability

Unique features such as the fixation targets, and the 1mm square light patch for assessing a/c flare set the new standard for Portable Slit Lamps.

#### For the digital age

Simply connect an iPhone 4 and the PSL becomes a portable, digital slit lamp. What could be simpler or more cost efficient?





#### Slit wheel

0	•	<b>O</b>	<b>O</b>	<b>O</b>	O
1x1	0.15	0.5	0.8	1.6	12

#### Filter wheel

Red free, Blue, Neutral density (0.8) and Clear

### Magnification

#### Slit and filter system Rotating Slit wheel selection Slit type Slit Lamp Slit length 12mm Slit widths 0.15mm, 0.5mm, 0.8mm and 1.6mm slits, 12mm circle and a 1mm square **Filters** Red free, Blue, Neutral density 0.8 and Clear IR protection Inbuilt IR cut filter $+/-60^{\circ}$ Slit angle Illumination control Continuously variable from low to full brightness **Fixation targets** LED system, double click of trigger **Power supply** Switch mode, (110V - 240V) Power supply unit +/- 10% multi-plug Power supply output 30VA (12V DC 2.5A) **Biomicroscope** Туре Binocular Hand held Biomicroscope Slit Lamp Optics Converging binoculars at 13° Magnification 10x and 16x, lever change Objective lens working distance @ 10x 16mm **Ordering details PSL Classic** 3010-P-2000 iPhone 4 adapter 3010-P-7010



3010-P-7000

Carrying case

Keeler Limited, Clewer Hill Road, Windsor, Berkshire SL4 4AA

+44 (0)1753 857177 | www.keeler.co.uk | info@keeler.co.uk

Keeler Instruments Inc. 456 Parkway Broomall, PA 19008 USA

Toll Free; 1 800 523 5620 | Fax: 1 610 353 7814 | www.keelerusa.com